

InsulScan®

Insulation Testing System

Prevent insulation failure at critical moments



When safety is a priority, trust InsulScan® to protect your patients, surgeons and hospital.

Specifically designed to test the insulation of electrosurgical instrumentation right in the operating room, InsulScan is safety meeting convenience. A quick check of the insulated instruments with InsulScan prior to the case will expose defective instruments so they can be removed and re-insulated. One last scan post-operatively will confirm that no defects occurred during surgery.

InsulScan Advantages:

- Designed for use in the Operating Room and Sterile Processing Department.
- Tests instrument pre-operatively and post-operatively.
- Available for sterile and non-sterile use.
- Checks both disposable and reusable instruments.

Safety meets Convenience

InsulScan is a reliable and convenient method for detecting failures in instrument insulation.

Each InsulScan system is comprised of:

- 9400-02 Insulscan Control Unit
- 9400-03 Sterile Wands (Box of 50)
- 9400-03U50 Non-Sterile Wands (Box of 50)
- 9400-05 Alligator Clips (Set of 2)
- 9400-06 Battery Charger
- 9400-07 Insulation Testing Brush
- 9400-08 Power Cord
- 9400-09 Reusable Cable

To order InsulScan contact your local representative or call 800.722.3675