

# Legacy/APM™ Control Unit

Alternation Therapy Air Pump



## Key Specs

### CONTROL UNIT

Cycle Times ..... 10, 15, 20, 25 min.  
Alarm ..... Low Pressure  
Pressure Sensor Technology ..... Yes  
Unit Weight ..... 4.5 lbs.

### MATTRESS OPTIONS

Legacy/APM™ Control Unit available with NPT3™, NP12™, and NP9™ non-powered mattresses.

The Legacy/APM Control Unit is an option for Sizewise NPT3™, NP12™, and NP9™ foam mattresses. The pump provides Static and Alternation therapy and is fully adjustable for patient comfort and healing. Simple controls include cycle times (10 to 25 minutes), soft/firm, and Autofirm max-inflate function. Control unit hangs from the footboard of most medical bed frames.

- Static and 1-in-2 Alternation therapy air pump
- 10, 15, 20, or 25-minute Alternation cycle times
- Whisper quiet with adjustable pressure settings for patient comfort
- Constant monitoring of air pressure; low pressure alarm and LED

[sizewise.com](http://sizewise.com) 800-814-9389

