

Bari Rehab Platform 2™

Bariatric Med-Surg/ICU Bed



Shown with standard pressure redistribution mattress and optional Power Drive

Bari Rehab Platform 2™ (BRP2) is a versatile 1,000 lb. capacity bed which navigates easily from Med-Surg to ICU to Rehab. BRP2 comes standard with an expandable pressure redistribution mattress, but its open architecture flat deck supports most surfaces—Sizewise or other brands. Numerous features help patients become actively involved in their care plan.

- *Simplified hand control and stand-assist siderails increase patient independence*
- *Head angle indicator assists with ventilator-associated pneumonia (VAP) protocols*
- *Optional Power Drive™ for easy navigation*
- *Available with Hercules Patient Repositioner™*
- *Low working deck height: 15"*

sizewise.com 800-814-9389


Sizewise®
Designed to heal



- 1 — Rest Secure System™ includes bed exit alarm, protocol timer, and scale
- 2 — Optional Power Drive™ for easy navigation
- 3 — Available with Hercules Patient Repositioner™—reposition with the touch of a button
- 4 — Expandable deck easily slides from 39"-48"

Key Specs

BED

| | |
|------------------------------|------------|
| Low Position: Deck to Floor | 15" |
| High Position: Deck to Floor | 30" |
| Deck Width | 39-48" |
| Deck Length | .80 or 86" |
| Safe Working Load | 1,000 lb. |

SIDERAILS

| | |
|----------|--------------------------------|
| Type | Removable 3-Position Fold Down |
| Included | 2 Standard; Foot Optional |

HEAD/FOOTBOARDS

| | |
|----------|--------|
| Material | Kydex® |
| Included | Yes |

CASTER AND BRAKING SYSTEM

| | |
|-------------------------|----|
| Dual Brake/Steer System | 5" |
|-------------------------|----|

BED ANGLES

| | |
|-----------------------|----------|
| Backrest/Head Section | 0 to 60° |
| Knee Section | 0 to 50° |
| Trend | 15° |
| Reverse Trend | -15° |

OPTIONS AND ACCESSORIES

| | |
|---|----------|
| Rest Secure System™ (RSS) | Standard |
| <i>Bed Exit Alarm, Protocol Timer, Built-in Scale</i> | |
| Battery Backup | Optional |
| IV Pole | Optional |
| O2 Holder | Optional |
| Pillow Speaker | Optional |
| Power Drive™ | Optional |
| Trapeze | Optional |
| 6-Port Aux Outlet | Optional |

MATTRESS OPTIONS

Standard pressure redistribution mattress included with purchase. Numerous non-powered and low air loss therapeutic surfaces available.

| | |
|-----------------------------------|----------|
| Hercules Patient Repositioner | Optional |
| Air Hercules Patient Repositioner | Optional |

size wise.com 800-814-9389

