

Bari Rehab Platform 3™

Bariatric Med-Surg/ICU Bed



Shown with optional pressure redistribution mattress

Bari Rehab Platform 3™ (BRP3) debuts the Power Drive 360™ for easy navigation through crowded hospital corridors—forward, backward, and side-to-side, plus the ability to turn on a dime. Numerous features help patients become more actively involved in their care plan.

- *Optional Power Drive 360™ for easy navigation*
- *Open architecture flat deck accepts most mattresses—Sizewise or other brands*
- *Fully compliant with latest IEC safety standards (60601-2-52)*
- *Available with Hercules Patient Repositioner™*

sizewise.com 800-814-9389


Sizewise®
Designed to heal



- 1 — Optional Power Drive 360™ moves forward, backward, and side-to-side with ease
- 2 — Rest Secure System 3™ includes bed exit alarm, protocol timer, and scale
- 3 — Open architecture flat deck accepts most mattresses—Sizewise or other brands
- 4 — Available with Hercules Patient Repositioner™

Key Specs

BED

Low Position – Deck to Floor	15.35"
High Position – Deck to Floor	30"
Deck Width	39-48"
Deck Length	80 or 86"
Safe Working Load	1,000 lb.
Expandable Deck	Yes

SIDERAILS

Type	Removable 2-Position Fold Down
Included	2 Standard; Foot Optional

HEAD/FOOTBOARDS

Type	Kydex®
Accessory Bar	Standard
Integrated Info Holder	Standard

CASTER AND BRAKING SYSTEM

Dual Brake/Steer System	5"
-------------------------	----

BED ANGLES

Backrest	0 to 64°
Knee Section	0 to 54°
Trend	15°
Reverse Trend	-15°

OPTIONS AND ACCESSORIES

Nurse Call Interface	Standard
Rest Secure System 3™ (RSS)	Standard
Bed Exit Alarm, Protocol Timer, Scale	
RSS Other Languages	Optional
Battery Backup	Optional
IV Pole	Optional
Line Management	Optional
O2 Holder	Optional
Pillow Speaker	Optional
Trapeze	Optional
6-Port Aux Outlet	Optional

MATTRESS OPTIONS

Numerous non-powered and low air loss therapeutic surfaces available.

Hercules Patient Repositioner	Optional
Air Hercules Patient Repositioner	Optional

sizewise.com 800-814-9389


Sizewise®
 Designed to heal