

# Behavioral Health Bed™ – Platform

Ligature-Resistant Bed



*Shown with optional low-profile head/footboard and lift opening*

Behavioral Health Bed™ – Platform provides optimal patient and caregiver safety. Made from high-density polyethylene and built with tamper-resistant fasteners, the Platform is durable and easy to clean. Anchor points fasten the bed to the floor, ensuring the bed is stable and secure. Available in two colors (Mist or Cocoa) with multiple head/footboard profiles, lift openings, and bed risers.

Featured in the  
**Behavioral Health  
Design Guide**<sup>1</sup>

- *High-density polyethylene: durable and easy to clean*
- *Tamper-resistant fasteners*
- *Anchor points allow bed to be secured to floor*
- *Safe working load of 1,000 lb.*
- *Full-enclosure base and low-profile head/footboards available*

size-wise.com 800-814-9389

  
**Size-wise**<sup>®</sup>  
Designed to heal



## Key Specs

**BED**

Sleep Surface Width	35"
Sleep Surface Length	82"
Platform Height	14"
Safe Working Load	1,000 lb.
Product Weight (open base)	250 lb.

*Note: Low-Profile Head/Footboards reduce bed weight by approximately 40 lb.*



**OPTIONS**

Anchor Points, Floor	Standard
Anchor Points, Wall*	Optional
Colors	Mist or Cocoa
Bed Riser Base	Optional
Custom Head/Footboard	Optional
Full-Enclosure Base	Optional
Lift Opening	Optional
Low-Profile Head/Footboard	Optional
Restraint Slots	Optional

\*Custom anchor points available upon request and approval by Sizewise.



**MATTRESSES**  
 Behavioral Health Mattress™ available with RF-welded seams for safety.

- 1 — Two color options: Mist (see front) or Cocoa (shown with standard head/footboard and optional restraint slots)
- 2 — Optional lift opening accommodates patient lifts
- 3 — Four anchor points allow bed to be secured to floor
- 4 — Available Behavioral Health Mattress™ designed with safety in mind
- 5 — Optional restraint slots

[sizewise.com](http://sizewise.com) 800-814-9389

