

Horizon 13.1™

Long-Term Care Bed Frame



Horizon 13.1™ is a long-term care bed frame with an expandable/retractable deck for accommodating most patient sizes. An adjustable wall saver protects walls from bed edges, and synchronized high-low keeps the bed level for positioning quietly and smoothly.

- Expandable/retractable deck width: 35" to 42"
- Adjustable wall saver to protect walls
- Synchronized high-low to keep bed level
- Easy-to-use 2-pedal caster locking system
- Optional trapeze kit

Key Specs

BED

Low Position – Deck to Floor	7"
High Position – Deck to Floor	30"
Deck Width	35"-42"
Deck Length	76" or 80"
Safe Working Load	600 lb.

SIDERAILS

Premium Counter-Rotating Assist Device	Optional
Premium Rotating Assist Bar	Optional
Premium Fixed Assist Bar	Optional

HEAD/FOOTBOARDS

Shape	Square
Included	Yes

BED ANGLES

Trend	Standard
Reverse Trend	Standard
Chair	Standard

OPTIONS AND ACCESSORIES

Adjustable Length (84"-88")	Optional
Enhanced Embedded Staff Control	Optional
Embedded Staff Control	Optional
Pan Deck Kit	Optional
Trapeze Kit	Optional
USB Power Supply	Optional

sizewise.com 800-814-9389

