

NP Duo™

Two-Layer Foam
Support Surface



NP Duo features two layers of foam to provide a consistent level of support and immersion from head to heel. The top layer of visco foam is softer for comfort, while the firmer bottom layer adds stability and prevents bottoming out. The surface also includes highly chemically resistant performance fabrics and RF-welded seams to improve durability and minimize risk of fluid ingress.



Two layers of foam provide universal head-to-heel pressure redistribution



4-way stretch, high-performance breathable top cover promotes immersion



RF-welded seams improve durability and minimize risk of fluid ingress

KEY SPECS

Max Patient Weight 500 lb.
Width 35"
Length 82"
Height 7"

CONSTRUCTION DETAILS

Foam Layers 2
RF-Welded Seams Standard

OPTIONS

Behavioral Health NP Duo BH™
Expandable NP Duo EX™
Custom Sizes By Request

LEARN MORE

Ask your Agiliti representative or call 800-814-9389