

NP Trio™

Three-Layer Advanced
Foam Support Surface



CuPro™ foam

NP Trio features three layers of foam in three zones to support varying pressure redistribution needs in the head, torso and heel sections. A four-way stretch top cover fabric and a layer of CuPro copper- and gel-infused foam help enhance microclimate by dissipating heat and moisture. The surface also includes highly chemically resistant performance fabrics, RF-welded seams and exclusive Core Shield™ layer to protect the internal components from contamination.



Three layers and three zones of foam for targeted pressure redistribution



4-way stretch, high-performance breathable top cover promotes immersion



Core Shield minimizes risk of fluid ingress, extending surface life

KEY SPECS

Max Patient Weight 500 lb.
Width 35"
Length 82"
Height 7"

CONSTRUCTION DETAILS

Foam Layers 3
Edge Support Standard
RF-Welded Seams Standard

OPTIONS

Behavioral Health NP Trio BH™
Seamless Expandable NP Trio EX™
Custom Sizes By Request

LEARN MORE

Ask your Agiliti representative or call 800-814-9389

