

# Adapt Pro<sup>Elite</sup>

Alternating Pressure Enhanced Low Air Loss Mattress System with Patient Sensing Technology

## MODEL 9200 - ADAPT PRO ELITE™ WITH PATIENT SENSING TECHNOLOGY

The Adapt Pro Elite™, Alternating Pressure Enhanced Low Air Loss Mattress with Patient Sensing Technology (PST), is the latest technology for the customized treatment and prevention of stages I-IV pressure sores. The Adapt Pro Elite™ Patient Sensing Technology (PST) measures the patient and automatically adjusts to create the patient's optimal therapeutic surface. Low Air Loss is achieved through a unique Ventilated Low Air Loss Cover that enhances therapy and helps to control micro-climate. Designed to fit all standard hospital beds, Adapt Pro Elite™ supports patient weights up to 400 lbs. The user friendly controls and ease of setup make Adapt Pro Elite™, ideal for institutional use and the home-bound patient. Adapt Pro Elite™ meets California Technical Bulletin #117 for fire retardancy.

CPR & Bolster Release valves quickly deflate mattress system for medical emergencies.



CPR



Bolster Release

Adapt Pro Elite™ with air safety bolster.



Shown: Model 9200SR5

### OPTIONS

Optional Electric Inflator



## DURABLE DUAL BLADDERS DESIGN

20 independent 8" deep air cells is constructed of (TPU) thermal polyurethane, laminated to a medical grade nylon for maximum durability and performance. The dual bladder design maintains air and prevents patients from bottoming out during power outages.

## VYVEX-I® ENHANCED LOW AIR LOSS COVER

Low Air Loss therapy is incorporated into the mattress cover, transferring moisture away from the patient's skin and maintaining the optimal micro-climate. Vyvex-I® is a low shear, low friction nylon cover fabric that is washable, anti-microbial, fluid resistant and fire retardant. Cover is loose fitting and will conform to the patient.

## PUMP CONTROL MODES & FEATURES

Adapt Pro Elite's dual compressor pump is feature rich, whisper quiet and extremely durable.

### MODES

- Alternating
- Static
- Static Care - falls back to Alternating in 20 mins.

### AUTOMATIC PATIENT SENSING TECHNOLOGY

- Auto Detection (PST)
- Fowler - Low 25° & Standard 45°

### FEATURES

- PST Patient Sensing Technology
- Comfort Control
- Cycle Time 10 - 30 minutes
- Auto-Firm
- Easy to read digital display
- Low Pressure & Power Failure Alarms

Blue Chip

Medical Products, Inc.

"Setting the Standard"

# ADAPT PRO ELITE™

## SPECIFICATIONS

### MATTRESS ONLY

Dimensions:	36" x 80" x 8"
Weight:	18 lbs.
Air Cells:	20, 8" TPU Nylon Air cells 3 Static cells in head section
Low Air Loss Cover:	Vyvetex-1® removable zippered, low shear, moisture resistant / vapor permeable cover

Adjustable anchor straps secure mattress to bed frame. Fire evacuation handles standard



### PUMP

Dimensions:	16.5" x 9" x 5"
Weight:	10 lbs.
Electrical:	120V 60Hz
Cycle Time:	10 - 30 minutes
Air Flow:	9 + 9 = 18 LPM
Listings:	UL CE CB
Ground Resistance:	0.1 ohms maximum
Current Usage:	less than 100 micro amps
Power Cord:	12 ft.
Construction:	ABS fire retardant UL 94V-0



## ADAPT PRO ELITE™ ORDERING INFORMATION

Model: 9200SR5	36" Mattress with Side Bolsters & Pump Control Unit
Model: 9200SR5-42	42" Mattress with Side Bolsters & Pump Control Unit
Model: 9200SR5-48	48" Mattress with Side Bolsters & Pump Control Unit
Model: 9200	36" Mattress without Side Bolsters & Pump Control Unit
Model: 9200-42	42" Mattress without Side Bolsters & Pump Control Unit
Model: 9200-48	48" Mattress without Side Bolsters & Pump Control Unit
Model: 9201	Pump Control Unit for 36" Mattress with or without Side Bolster
Model: 9201-42	Pump Control Unit for 42" Mattress with or without Side Bolster
Model: 9201-48	Pump Control Unit for 48" Mattress with or without Side Bolster
Model: 9220SR5	36" Mattress only with Side Bolsters
Model: 9220SR5-42	42" Mattress only with Side Bolsters
Model: 9220SR5-48	48" Mattress only with Side Bolsters
Model: 9202	36" Mattress only without Side Bolsters
Model: 9202-42	42" Mattress only without Side Bolsters
Model: 9202-48	48" Mattress only without Side Bolsters
Model: 9623	Electric Inflator & Deflator

Adapt Pro Elite™ meets Group 2 Support Surfaces requirements for Medicare E0277. Patient must meet all medical criteria to qualify for reimbursement.

### WARRANTY

Mattress & Pump Cover:	1 year limited 6 months limited
------------------------	------------------------------------



7-11 Suffern Place, Suffern NY 10901

Tel: (845) 369-7535 • Fax: (845) 369-7633 • Toll Free: (800) 795-6115

www.bluechipmedical.com

"Setting the Standard"

©Blue Chip Medical Products, Inc. 2019

Alternating Pressure Mattresses  
Alternating Pressure Pads & Pumps  
Low Air-Loss Mattress Systems  
Lateral Rotation Mattress Systems  
Gel Mattresses  
Gel Mattress Overlays  
Static-Air Mattress Overlays  
Water Mattress Overlays  
Wheelchair Cushions  
Patient Lift Slings  
Side Rail Bumpers  
Cart Covers  
OR Replacement Table Pads  
Replacement Stretcher Pads  
Cubicle Curtains  
Hospital Pillows