

# Bedside Recliner™

Comfort and Mobility Item



Bedside Recliner is designed for both comfort and safety. Features include three swivel-locking casters and one rigid-locking caster, plus a moisture barrier over the high-density foam that allows for easy cleaning. Calf and footrest combined with multiple reclining positions provide optimal patient comfort.

- *High-density foam covered in moisture barrier for easy cleaning*
- *Swivel- and rigid-locking casters for flexibility and safety*
- *Breakaway footrest and ergonomic push bar*
- *Quick-release feature allows for removable back*
- *Standard and bariatric versions available*

size wise.com 800-814-9389

  
**Size wise®**  
Designed to heal



## Key Specs

Weight Capacity (Standard) .....	500 lb.
Weight Capacity (Bariatric) .....	750 lb.
Seat Width (Standard) .....	24"
Seat Width (Bariatric) .....	29"
Seat Depth .....	21"
Seat Height .....	21.5"



- 1 — Features swivel- and rigid-locking casters for flexibility and safety
- 2 — Steel-to-steel construction for longer life
- 3 — Ergonomic push bar

size wise .com 800-814-9389

