

# Behavioral Health Bed™

Fully Electric, Ligature-Resistant Bed



The first of its kind, Behavioral Health Bed™ is a fully electric bed frame designed in collaboration with behavioral health clinicians to meet the needs of dual-diagnosis patients. The bed ensures patient safety by using PowerSafe™ tamper-resistant fasteners on all electrical parts and lockout controls accessible with a universal key. An anchor point allows the bed to be secured to the floor.

Featured in the  
**Behavioral Health  
Design Guide¹**

- Fully electric with easy-to-use footboard lockout controls
- Electrical articulation for head, knee, and high/low systems
- Contiguous molded head/footboards and siderails
- Tamper-resistant features include encased wiring and security screws
- Designed in collaboration with behavioral health clinicians

[sizewise.com](https://sizewise.com) 800-814-9389

  
**Sizewise®**  
Designed to heal



- 1 — Tamper-resistant fasteners and security screws
- 2 — Keyed footboard lockout for added safety
- 3 — Optional locking storage box stows removable power cord
- 4 — Optional Rest Secure System™ with bed exit alarm, protocol timer, and scale

## Key Specs

### BED

|                                      |            |
|--------------------------------------|------------|
| Low Position - Deck to Floor         | 7.5"       |
| Low Position - Deck to Floor (w/RSS) | 9"         |
| High Position - Deck to Floor        | 30"        |
| Deck Width                           | 35" or 39" |
| Deck Length                          | 82"        |
| Safe Working Load                    | 850 lb.    |
| Bariatric Version                    | Available  |

### SIDERAILS

|              |            |
|--------------|------------|
| Type         | Swing Down |
| Head Section | Standard   |
| Foot Section | Optional   |

### HEAD/FOOTBOARDS

|      |                     |
|------|---------------------|
| Type | High-Impact Plastic |
|------|---------------------|

### CASTER AND BRAKING SYSTEM

|                         |    |
|-------------------------|----|
| Caster Diameter         | 4" |
| Dual Brake/Steer System | 4" |

### BED ANGLES

|               |          |
|---------------|----------|
| Head Section  | 0 to 58° |
| Knee Section  | 0 to 49° |
| Foot Section  | 0 to 8°  |
| Trend         | 15°      |
| Reverse Trend | -15°     |

### OPTIONS AND ACCESSORIES

|                                       |          |
|---------------------------------------|----------|
| Battery Backup                        | Standard |
| Keyed Footboard Lockout Controls      | Standard |
| Quick-Disconnect Hand Control         | Standard |
| Removable Power Cord                  | Standard |
| Tamper-Resistant Fasteners            | Standard |
| Behavioral Health Mattress™           | Optional |
| Deck Gap Fillers                      | Optional |
| Floor Safety Mat                      | Optional |
| Locking Storage Box Assembly          | Optional |
| Pedal Lock Kit                        | Optional |
| Restraint Slots                       | Optional |
| Rest Secure System™ (RSS)             | Optional |
| Bed Exit Alarm, Protocol Timer, Scale |          |
| Seizure Pad Covers                    | Optional |
| Color: Mist                           | Standard |
| Color: Cocoa                          | Optional |

sizewise.com 800-814-9389

