

# Carewise 9.1™

## Foam Support Surface



### Key Specs

Capacity .....	500 lbs.
Width .....	35", 39", or 42"
Length .....	76", 80", or 84"
Height.....	7"

### CONSTRUCTION

Foam Layers.....	5
Defined Perimeter.....	Optional

### FABRIC OPTIONS

4-Way Stretch Top Cover .....	Standard
PVC Laminate Bottom Cover .....	Standard
Fire Barrier .....	Standard

Carewise 9.1™ features seven pressure redistribution zones housed in viscoelastic memory-style foam for optimal patient comfort. A diamond-cut top foam layer for pressure redistribution, shear reduction, and air circulation, keeping the surface cool. The heel area is softer than the rest of the surface to protect vulnerable heels, and a concealed zipper prevents skin tears and protects inner foam.

- Seven pressure redistribution zones
- Diamond-cut top foam layer for cool surface
- Heel area softer than rest of surface
- Concealed zipper prevents skin tears and protects inner foam
- Made-to-order item

size wise.com 800-814-9389

