

Reliable Solutions for Critical Access Hospitals

Critical Access Hospitals need beds and support surfaces that are flexible enough to serve multiple patient acuities, durable enough to ensure high uptime and utilization, and budget-friendly enough to keep service and repair bills to a minimum.

Sizewise offers a one-size-fits-most combination of Alliance® Med-Surg bed frame and NPT3™ pressure redistribution support surface. The products are easy to use, require little maintenance, and have a track record of reliable service in hospitals of all sizes.



Alliance®
Med-Surg Bed



NPT3™
Pressure Redistribution Support Surface

Key Features/Benefits	<ul style="list-style-type: none"> • Acuity-adaptable to CAH and Swing Bed needs: Med-Surg to Rehab/Post-Acute • Ultra-lightweight and highly maneuverable • Open-architecture flat deck accepts most mattresses: Sizewise and other 	<ul style="list-style-type: none"> • Advanced medical fabric stands up to hospital cleaning protocols • 12 foam-filled cells provide head-to-heel support for optimal patient comfort • Optional 1-in-2 alternation therapy adds higher level of pressure injury prevention
Warranty*	Frame: 5-Year Limited Electronics and Parts: 2-Year Limited	Construction: 5-Year Limited Fabric and Foam: 2-Year Limited Electronics/MAP Liner: 1-Year Limited

**See product user manual for a full copy of warranty terms and conditions.*

The Sizewise Promise

- Budget-friendly capital, lease-to-own, and rental options
- Simplicity of design: easier to use, less likely to need costly repairs
- 24/7/365 customer service and support

Learn more: Call 800-814-9389
or contact your Sizewise rep

