

Electri-Cool® II

Localized Cold Therapy System



THE ELECTRI-COOL® II IS THE PREMIER CHOICE OF LOCALIZED COLD THERAPY SYSTEMS.

- Uses NO ICE, which means no messy wet floor or worry about using valuable nursing time to change the ice.
- Helps decrease pain and could reduce the amount of medication needed for the patient
- More consistent cooling therapy

Possible applications for the Electri-Cool® II range from surgical incision treatment, swelling from injury, sprains, strains, and more. The Electri-Cool® II is lightweight and compact, making this unit the ideal choice among healthcare professionals. It may also be attached to the patient's bed with a bracket to save valuable floor space or may be mounted on a stand for easy mobility.

Cincinnati Sub-Zero has been managing your patient's temperature needs for over 45 years.

CSZ
Cincinnati Sub-Zero

Electri-Cool® II

Localized Cold Therapy System

| Ordering Information | | |
|------------------------------------|---|-------------------------------|
| Cat. No. | Description | Quantity |
| 767 | Electri-Cool® II | 100v/115v, 230v, 220v/240v |
| 767-BBK | Bed Bracket for Electri-Cool® II | ea. |
| UNV LPS | Universal Stand | ea. |
| EC-080 | Adapter Plate for Electri-Cool® II (for use with universal stand) | ea. |
| TF-100 | Dual-Matic Flow Rate and Temperature Tester | ea. |
| Temp-Pads® Cold Therapy Pad | | |
| 510-US | Temp-Pad® - 5" x 10' (12.7cm x 25.4) Sterile | 10/carton |
| 510-USS | Temp-Pad® - 5" x 10" (12.7cm x 25.4cm) Sterile with straps | 10/carton |
| 814-US | Temp-Pad® - 8" x 14" (20.3cm x 35.6cm) Sterile | 10/carton |
| 814-USS | Temp-Pad® - 8" x 14" (20.3cm x 35.6cm) Sterile with straps | 10/carton |
| 1112-US | Temp-Pad® - 11' x 12" (27.9cm x 30.5cm) Sterile | 10/carton |
| 1112-USS | Temp-Pad® - 11" x 12" (27.9cm x 30.5cm) Sterile with straps | 10/carton |
| 1215-US | Temp-Pad® - 12" x 15" (30.5cm x 38.1cm) Sterile | 10/carton |
| CT-99 | Temp-Pad® - Special Knee/Shoulder pad | 10/carton |
| ImmobilICE Supports (includes pad) | | |
| 814-IMB | Back support | 5/carton |



Product Specification:

Dimensions:

12" W x 6" D x 12" H
(30.5cm x 15.2cm x 30.5cm)

Weight: 13 lbs. (5.9kg)

Water Temperature Range:

40°F-45°F (5°C-7°C)
45°F-50°F (8°C-10°C)
50°F-55°F (11°C-13°C)

CSZ
Cincinnati Sub-Zero

Phone: 513-772-8810
Toll Free: 800-989-7373
Fax: 513-772-9119

www.cszmedical.com



© 2009 by Cincinnati Sub-Zero Products, Inc.
99194 4/10