

Envy® S Series™

Pressure Redistribution Stretcher Pad





The Envy® S Series™ stretcher pad system features F3 Free™ construction, which reduces harmful chemicals in the fabric, foam, and fire barrier. The pad contours to the patient's body and provides low-interface pressure, aiding in pressure injury prevention. A high-density foam perimeter provides additional stability while enhancing comfort.

- F3 Free™ construction reduces harmful chemicals in the fabric, foam, and fire barrier
- Provides low-interface pressure and optimal contour
- High-density foam perimeter for additional stability
- Heel slope offloads pressure from the heels for extra comfort

size wise.com 800-814-9389


Size wise®
Designed to heal



1 — Pressure Mapping: 50 mmHg scale shown

2 — Heel slope for optimal pressure relief

Key Specs

PAD

Capacity 500 lbs.
 Width Custom
 Length Custom
 Solid Foam Mattress Available

FABRIC OPTIONS

Envy® 4-Way Stretch Top Cover Standard
 Extreme Polycarb 4-Way Stretch Top Optional
 Envy® Non-Skid Bottom Cover Standard
 Halogen-Free Fire Barrier Standard
 MAP Liner™ Optional

FABRIC CHARACTERISTICS



4-Way
Stretch



Chemical
Resistant



Green



A PRACTICE GREENHEALTH PROGRAM

Sizewise proudly supports Practice Greenhealth's Healthier Hospitals Safer Chemicals Challenge for Healthy Interiors. Our Envy products meet or exceed Healthier Hospitals goals of reducing patient and caregiver exposure to chemicals of concern.

sizewise.com 800-814-9389



Sizewise®
 Designed to heal