

Evolution[®]

Fall Prevention Bed



Shown with optional defined perimeter foam support surface

Evolution is a go-anywhere fall prevention bed with a 9" deck low and 30.75" deck high for an ergonomic work environment. The bed transports patients in any high or low position. Rest Secure System[™] (RSS) comes standard. Combine Evolution with a low air loss surface to provide alternation, pulsation, rotation, or immersion therapy.

- Focus on fall prevention: 9" deck low
- Rest Secure System[™] (RSS) includes bed exit alarm, protocol timer, and scale
- Transports in any high or low position
- Evolution Expandable[™] accommodates 42" or 48" mattresses

size wise.com 800-814-9389


Size wise[®]
Designed to heal



Key Specs

BED

Low Position: Deck	.9"
High Position: Deck	30.75"
Deck Width (Standard)	35" or 39"
Deck Width (Expandable)	39", 42", 48"
Deck Length	82"
Safe Working Load	.850 lbs.

SIDERAILS

Type	Fold Down
Included	2 Standard; Foot Optional

HEAD/FOOTBOARDS

Type	High-Impact Plastic
------	---------------------

CASTER AND BRAKING SYSTEM

Caster Diameter	.4"
Dual Brake/Steer System	.4"

BED ANGLES

Head Section	.0 to 58°
Knee Section	.0 to 49°
Trend	15°
Reverse Trend	-15°

OPTIONS AND ACCESSORIES

Battery Backup	Standard
CPR Release	Standard
Floor Safety Mat	Standard
Footboard Lockout Controls	Standard
Rest Secure System™ (RSS)	Standard
<i>Bed Exit Alarm, Protocol Timer, Built-in Scale</i>	
IV Pole	Optional
O2 Holder	Optional
Pillow Speaker	Optional
Seizure Pad Covers	Optional
Trapeze	Optional
6-Port Aux Outlet	Optional
Color: Mist	Standard
Color: Cocoa	Optional

MATTRESS OPTIONS

Numerous non-powered and low air loss therapeutic surfaces available.

1 — The standard Evolution® has a deck low of 9", the lowest full-feature bed frame. Also available, Low Boy®, with an industry-leading 7.5" low deck

2 — Optional expandable deck: adjusts from 39" to 42" or 48"

3 — Rest Secure System™ includes bed exit alarm, protocol timer, and scale

size wise.com 800-814-9389

