

Low Boy®

Fall Prevention Bed




Industry-Leading 7.5" Deck Low

Low Boy® fall prevention bed features the *lowest* deck low on the market at 7.5", plus the ability to transport in any high or low position. The open architecture flat deck accepts most mattresses, and the bed features lockout controls for patient safety and caregiver peace of mind. Floor Safety Mat included with purchase.

- *Focus on fall prevention: 7.5" deck low*
- *Bariatric and Expandable Deck versions are available*
- *Transports in any high/low position*
- *Open architecture flat deck accommodates most mattresses*
- *Lockout controls for patient safety and caregiver peace of mind*

sizewise.com 800-814-9389


Sizewise®
Designed to heal



- 1 — Open architecture flat deck accommodates most mattresses
- 2 — Lockout controls on footboard
- 3 — Expandable deck option allows for 42"- or 48"-wide support surface

Key Specs

BED

Low Position – Deck to Floor	7.5"
High Position – Deck to Floor	30.75"
Deck Width (Standard)	35" or 39"
Deck Width (Expanded)	39", 42", 48"
Deck Length	82"
Safe Working Load	850 lbs.

SIDE RAILS

Type	Fold Down
Included	1 Set

HEAD/FOOTBOARDS

Type	High-Impact Plastic
Included	Yes

CASTER AND BRAKING SYSTEM

Caster Diameter	4"
Dual Brake/Steer System	5"

BED ANGLES

Head Section	0-64°
Knee Section	0-54°
Foot Section	0-10°
Trend	15°
Reverse Trend	-15°

OPTIONS AND ACCESSORIES

Battery Backup	Standard
Floor Safety Mat	Standard
Bed Alarm Kit	Optional
CPR Release	Optional
Expandable Deck	Optional
IV Pole	Optional
O2 Holder	Optional
Seizure Pad Covers	Optional
Trapeze	Optional
6-Port Aux Outlet	Optional

MATTRESS OPTIONS

Numerous non-powered and low air loss therapeutic surfaces available.

sizewise.com 800-814-9389

