


## Matteo Motorized Recliner

Accommodates patients up to 750 lbs.



Dimensions:  
Width:36" Depth: 39"  
Extended depth: 72.5"  
Seat Width: 29.5"  
Seat Depth: 22"

The Matteo motorized bariatric recliner features a motorized mechanism with infinite reclined positions and a lift feature. Its load capacity of 750 lbs is made possible by two heavy duty motors, and a sturdy steel inner frame. A very simple control wand with only two buttons makes this a very user friendly chair, easily operated by patients and staff.



## Matteo Multilock Recliner

Accommodates patients up to 750 lbs.



Dimensions:  
Width:36" Depth: 32"  
Extended depth: 67"  
Seat Width: 28"  
Seat Depth: 21"

With a load capacity of 750 lbs and a seat width of 28", the Matteo bariatric recliner possesses all the technical features required by bariatric patients, within a simple and discreet design that provides comfort and dignity to the user. It's independent foot rest mechanism is operated by a lever. The back is released using a finger pull and locks in any position from upright to fully reclined.

