

Mobilizing a Bariatric Patient

Sizewise Bariatric Solutions provide a level of support and comfort for the most challenging patients. Our products promote early mobility, patient dignity, and caregiver safety. From boosting patients in bed to transfer and transport, you'll find products that are Designed to heal™.




Lateral Transfer:

- Bari Rehab Platform™
- Air Transfer Device™
- Shuttle B Series Advanced™



Total Lift:

- Bari Lift and Transfer™
- Shuttle B Series Advanced™



Ambulate:

- Shuttle B Series Advanced™
- Bari Walker™
- Bari Chair™



Boost:

- Hercules Patient Repositioner™

Learn more: Call 800-814-9389
or contact your Sizewise rep

Bariatric Solutions



Bari Rehab Platform 2™ w/Hercules Patient Repositioner™

- 39" or 48"W x 80" or 86"L; 1,000 lb. capacity
- 15" deck low, 30" deck high



Bari Lift and Transfer™

- 750 or 1,000 lbs. capacity
- Automatic base width adjustability
- Battery powered with internal charger
- Scale optional



Bari Chair™

- 750 lbs. capacity
- 24" to 36" widths available
- Extended surface/capacity, elevating leg rests



Bari Walker™

- 750 lbs. capacity
- 27"W between supports
- Adjustable height: 31" to 40"



Shuttle B Series Advanced™

- 1,000 lbs. capacity
- Transforms from chair to stretcher
- Independently elevating leg rests
- Power tilt and recline for patient comfort and caregiver safety
- Stand assist poles



Air Transfer Device™

- 1,000 lbs. capacity
- Ideal for patient transfer
- Radiolucent for X-ray/MRI



Bari Drop-Arm Commode™

- 600/1,000 lbs. capacity



Bedside Recliner™

- 750 lbs. capacity
- 30"W x 21"D

Learn more: Call 800-814-9389
or contact your Sizewise rep