



umano medical

## AN ALL-IN-ONE HOSPITAL BED FOR MULTIPLE CARE ENVIRONMENTS

med-surg | bariatric care | palliative care



 **ook snow** » ALL + **Pilot Drive**™

One-step Infinite™ width expansion system, increased weight capacity, and low height all come together to improve staff efficiency, patient comfort, and safety.

Featuring Pilot Drive™ – the first powered drive system to offer an unmatched combination of low height, underbed clearance, and maneuverability.

### VERSATILITY • LOW HEIGHT • MOBILITY



**MEDICAL-SURGICAL**  
Low height, ease of use, and ease of cleaning makes this bed ideal for med-surg environments.



**BARIATRIC CARE**  
1060 lb weight capacity, width expansion, and an appealing design encourages patient autonomy and dignity for bariatric care.



**PALLIATIVE CARE**  
Also known as the Cuddle Bed, it is the only bed on the market certified to be used by two people.



**PILOT DRIVE**  
Added mobility combined with enhanced versatility makes the ook snow ALL with Pilot Drive and integrated O2 bottle holders a unique solution for your bariatric patients in critical care.



STANDARD FEATURES

- 11" low height with 5" dual casters
- Infinite™ width adjustment
- Integrated scale system
- Smart screen, nurse controls and patient control lockout on footboard
- Integrated bed exit detection system
  - 3 levels of sensitivity
  - Bed exit side view lights
- Smart Nightlights™
- Remote bed exit alarm connectivity
- Boostless™ Backrest System
- One hand release siderails
- Lifiable and flat deck covers
- Battery back-up
- Compatibility with mattresses within 48" width
- Programmable lowest height
- Nurse and patient controls integrated into siderails
- Integrated bed extender (35" x 4")
- CPR release handle
- Trendelenburg angle displayed on footboard
- Backrest angle indicator
- Patient restraint loops
- Integrated pump holder
- 120-volt / 5A insulated outlets (at foot end)
- Communication port (for bed software updates)
- Programmable Auto-Contour mode
- One-touch button groups
  - 30-degree head of bed
  - One-button chair position
  - Flat bed position
  - Trendelenburg / Reverse Trendelenburg
  - Vascular foot position

OPTIONAL FEATURES

- 12.7" low height 6" dual casters
- Pilot Drive™ with or without O2 bottle holders (compatible with size E bottles)
- Downward motion interrupt
- Bed status monitoring
- Optional patient surface length (80" or 84")
- 120-volt/5A insulated outlets with USB mobile device power outlets
- Higher siderails (for 8" mattress)
- Foot end back-up control on the frame
- Roller type wall protection bumpers
- Additional nurse call port
- One-touch button groups
  - Backlit patient controls
  - Nurse call
- Different styling options available

ACCESSORIES

- Patient helper system (300 lb capacity)
- Upright oxygen bottle holder
- IV poles (permanent & removable)
- 4" bed extender cushion
- Removable patient controls
- Transport cart

SPECIFICATIONS

ooksnow » ALL

Product model	200-0000
Length	Overall length (80" deck length) 90" / 228.6 cm Overall length (84" deck length) 94" / 239 cm Extended length Adds 4" / 10.1 cm with Pilot Drive Adds 1" / 2.5 cm with Pilot Drive / oxygen bottle holders Adds 2.25" / 5.7 cm
Overall width	Siderails down - 35" surface 39.8" / 101 cm Siderails up - 35" surface 41.5" / 105.4 cm Siderails down - 47" surface 51.5" / 130.8 cm Siderails up - 47" surface 52.8" / 134 cm
Safe Working Load	1060 lb / 481 kg
Height range (to deck top)	Low (with 5" / 12.7 cm dual casters) 11" / 27.9 cm High (with 5" / 12.7 cm dual casters) 30" / 76.2 cm Low (with 6" / 15.2 cm dual casters) 12.7" / 32.3 cm High (with 6" / 15.2 cm dual casters) 31.7" / 80.5 cm
Patient deck positioning	Backrest 0° to 60° Knee gatch 0° to 30° Trend. / Reverse Trend. ± 16° Boostless™ 5" / 12,7 cm
Patient deck dimensions	Width 35" to 47" / 89 cm to 119.4 cm Length 80" or 84" / 203 cm or 213 cm Extended Length Adds 4" / 10.1 cm
Caster diameter	5" or 6" / 12.7 cm or 15.2 cm with Pilot Drive 6" / 15.2 cm

Compliant with standards IEC 60601-1 and IEC 60601-2-52.  
Certified for use by two people (referring to one patient and one bystander.)

ELECTRONICS

Volts	120 V AC
Ampere rating	4.25A (9.25A with 120V outlets)
Current leakage	Compliant with IEC 60601-1
Frequency	60 Hz
Hospital grade plug	Yes
Battery type	Sealed rechargeable lead acid 1.2Ah, 2 * 12V
Battery type (Pilot Drive)	Sealed rechargeable lead acid 16.5Ah, 2 * 12V
Scale system accuracy	1% with a maximum difference of 2.2 lb / 1 kg (in all bed positions)

Umano Medical reserves the right to change specifications without notice.

Umano Medical Inc.

230, boulevard Nilus-Leclerc  
L'Islet QC G0R 2C0 CANADA  
Toll free number : 1 844-409-4030  
information@umanomedical.com  
umanomedical.com