

MEETING ALL CARE UNIT REQUIREMENTS

ARE YOU TIRED OF WAITING AROUND FOR A SPECIALITY BED?

The ook snow ALL bed makes admitting a patient a lot easier and provides dignity to the patient. Through its all-purpose design, hospital staff can quickly adapt the bed to accommodate a new patient without having to leave the room.

With its one-step width expansion system, increased weight capacity, and low height, the ook snow ALL offers game-changing standardization possibilities throughout all care settings, including bariatric, med-surgical, and palliative care units. Through a versatile design, you can now save on equipment handling time, clear your storage space, and eliminate rental costs.

Getting access to the right bed has never been so simple.

**CAREGIVER TESTIMONIAL**  
 "We no longer have to accustom ourselves to various specialty beds - it is now the bed that adapts to each of our patients!"



INFINITE™ WIDTH IN A SLEEK AND CLEAN LOOK

NOT YOUR AVERAGE BARIATRIC BED

The need to switch out to a specialty bed for bariatric patients is no longer necessary when you have the ook snow ALL bed.

With features such as Infinite™ width technology and integrated bed extender, staff can quickly adapt the beds to accommodate any patient size, providing them with the space and comfort they deserve.

Designed to appear like our med-surg bed, the ook snow ALL bed promotes patient satisfaction, dignity, and comfort, without compromising on low height and functionality.

**A BED THAT DOES IT ALL**  
 VERSATILE, EXPANDABLE, AND MOBILE!

The ook snow ALL with Pilot Drive™ is a unique powered drive system designed with mobility, versatility, patient safety, and product performance all-in-one.

Reimagined with ease of use and fall prevention in mind, the ook snow ALL with Pilot Drive maintains its exceptional low height while still providing essential underbed clearance.

Featuring intuitive visual prompts, the Pilot Drive system is easy to use.



BENEFIT FROM ALL THE FEATURES OF THE ook snow ALL



**LOW HEIGHT**  
 Exceptional low height reduces the risk of falls and facilitates patient mobility.



**WEIGHT CAPACITY**  
 70 lb to 1060 lb weight capacity makes the bed versatile for any patient size.



**INFINITE™ WIDTH ADJUSTMENT**  
 Patient deck can be expanded to accommodate support surfaces from 35" to 48".



**SAFE GAP MONITORING**  
 Software-integrated system monitors bed width to ensure safe gaps are maintained.



**HIGHLY EFFICIENT SCALE**  
 The integrated scale provides accurate weight in any bed or patient position.



**EASY TO CLEAN**  
 Open architecture, seamless smooth surfaces, and liftable deck covers facilitate cleaning procedures.



**BOOSTLESS™ BACKREST SYSTEM**  
 Backrest moves backward to reduce shearing and patient repositioning needs.



**FALL RISK MANAGEMENT**  
 Fall prevention protocols are strengthened with the combination of low height, bed exit technology, and ergonomic siderail design.



**BATTERY AUTONOMY**  
 Staff will feel confident covering long distances, even up ramps and steep slopes, thanks to a strong battery autonomy and powered drive release system.



**CONTROLLED DRIVING SPEED**  
 Variable speed capabilities will give staff full control over the desired speed and velocity, while helping to relieve excessive push and pull stress.



**MAXIMUM MANEUVRABILITY**  
 Smooth steering and sharp 360° turns are made effortlessly with the centralized 5<sup>th</sup> wheel.




**INTEGRATED OXYGEN BOTTLE HOLDER**  
 Safely transporting bariatric patients in intensive care units is made easy with optional integrated oxygen bottle holder.



**POWERED WIDTH EXPANSION SYSTEM**  
 Allows for width adjustment in one simple step.



**SMART NIGHTLIGHTS™**  
 Provides patients with a clear view of surroundings.



**LIFTABLE DECK COVERS**  
 Simplifies cleaning and maintenance procedures.



**POWER OUTLETS WITH USB PORTS**  
 Keeps patient mobile devices within reach to limit fall risks.



**SMART SCREEN**  
 Optimizes patient safety through advanced bed monitoring capabilities.





ALLOWING CLOSENESS BETWEEN LOVED ONES

A BED FOR TWO

Getting close to loved ones in care settings has not always been possible.

In palliative care environments, this much-needed proximity between patients and loved ones is now achievable with the ook snow ALL bed.

Often referred to as the "cuddle bed", the ook snow ALL is the only bed on the market certified for use by two people.



BELOVED TESTIMONIAL

"There are no words to describe that little bit of happiness and comfort for the time I got to spend with my husband. [...] This bed allowed me to be close to him without squatting like in the regular hospital beds."

- Brianna Miles\*, NL, Canada



AN ALL-IN-ONE HOSPITAL BED FOR MULTIPLE CARE ENVIRONMENTS

med-surgical | bariatric care | palliative care



LOW HEIGHT, UNDERBED CLEARANCE, AND MANOEUVRABILITY

A NEW SHIFT IN POWERED DRIVE SYSTEMS

The ook snow ALL with Pilot Drive™ is the first system to deliver on low height, full underbed clearance, and maximum manoeuvrability.

With its controlled driving force, variable speed capabilities, and intuitive controls, healthcare teams will feel confident and in control when moving patients during intrahospital transportation.



Visit our Virtual Showroom for an interactive 3D experience of our full line of hospital beds and accessories

Toll free number : 1 844 409-4030  
umanomedical.com

200-9079-ENG R1 Copyright Umano Medical 2022 - Printed in Canada

