

Incubator and Bassinet Pads

Foam Support Surface



Key Specs

MATTRESS

Capacity	100 lbs.
Width	Custom
Length	Custom

FABRIC OPTIONS

Polycarbonate 4-Way Stretch	Standard
4-Way Stretch (in Royal or Green)	Optional
Non-Skid Polyurethane Bottom Cover ...	Standard

Our Incubator and Bassinet Pads are latex-free and radiolucent for X-rays. The First Point of Contact™ top cover is free from chemicals of concern and is vapor-permeable, low-friction, and anti-shear. Bassinet pads are available with a rigid base. Available in standard and custom sizes.

- *Latex-free and radiolucent*
- *First Point of Contact™ top cover free of chemicals of concern, vapor-permeable, low-friction, and anti-shear*
- *Available with rigid base*
- *Standard and custom sizes to fit any incubator or bassinet*

size-wise.com 800-814-9389

PRD-0030C-0618

