

Pediatric Support Surfaces

Sizewise Solutions



Pediatric Surfaces


Pediatric patients suffer from pressure injuries (PIs) and skin breakdown issues, just like older patients, but products for adults are not designed to accommodate children. Babies and toddlers often sink too far into adult-sized low air loss surfaces. Finding pressure redistribution surfaces specifically for little ones is often a difficult task—*until now*.

Size-wise pediatric support surfaces are safe, comfortable, and—most importantly—designed to heal your smallest patients. Our top covers are latex-free, vapor-permeable, low-friction, and anti-shear, providing babies and toddlers the best PI prevention while giving nurses and parents peace of mind.


Pressure Mapping

Infants and children at increased risk of PI benefit from the right support surface. Pressure mapping shows how the right surface can lower interface pressures, helping prevent PIs. *All illustrations feature a 3-year-old male patient, 34.8 lb.*



Pediatric Pulse™




Crib Mattress with Gel



Crib Mattress with Memory Foam






Pediatric Pulse™


Crib Mattress

Incubator and Bassinet Pads


Weight Capacity	1.6-100 lb./35-200 lb.	100 lb.	100 lb.
Size(s)	Custom	Custom	Custom
Therapy Modes	Pulsation, Static	Reactive Pressure Redistribution	Reactive Pressure Redistribution
Key Features	<ul style="list-style-type: none"> Vapor-permeable, low-friction, and anti-shear Vented air chambers provide clinically effective low air loss Soothing comfort ideal for nonpharmacological treatment of Neonatal Abstinence Syndrome (NAS) Available in crib and Med-Surg sizes to fit a variety of patient needs 	<ul style="list-style-type: none"> Latex-free Radiolucent Halogen-free fire barrier First Point of Contact™ top cover, free of harmful chemicals and vapor-permeable, low-friction, and anti-shear 	<ul style="list-style-type: none"> Latex-free Radiolucent First Point of Contact™ RF-welded top cover: free of harmful chemicals and vapor-permeable, low-friction, and anti-shear Rigid base available
Rental Available	✓		



Features kid-friendly blower keypad design



Optional gel or memory foam topper



Custom sizes to fit any incubator or bassinet

Learn More
800-814-9389
sizewise.com/pediatrics

