

# Rest Secure System™

Multifunctional Patient Care Monitor



Available in green or blue  
(depending on bed frame)

Rest Secure System (RSS) incorporates a bed exit alarm, protocol timer, and integrated scale—all in one device. RSS is customizable to meet caregiver and patient needs, allowing for efficiency, safety, and comfort. Available on numerous bed frames—rental and capital.

#### — Bed Exit Alarm

- Three sensitivity settings; five alarm tones
- Voice-recordable alarm feature
- Interfaces with existing nurse call systems

#### — Protocol Timer

- Timer intervals: 1, 2, or 4 hours
- Audible tone and nurse call system interface

#### — Integrated Scale

- Weighs patient in bed; lbs. or kg
- Automatically weighs daily and stores 30-day history

## Key Specs

### RSS MODULE

Power . . . . . AC and Rechargeable  
Operating Time (Battery) . . . . . 200 min  
Recharge Time (Battery) . . . . . 30-120 min  
Voice Recording . . . . . up to 15 sec

### FUNCTIONS

Bed Exit Alarm . . . . . Off, Low, Med, High  
Protocol Timer . . . . . Off, 1, 2, 4 hrs  
Volume . . . . . Off, Low, Med, High  
Scale . . . . . 0-1000 lbs.

### PRODUCT AVAILABILITY

Bari Rehab Platform 2™ . . . . . Standard  
Bari Rehab Platform 3™ . . . . . Standard  
Evolution® . . . . . Standard  
Behavioral Health Bed™ . . . . . Optional

sizewise.com 800-814-9389

