

Sapphire Pulse™

Low Air Loss Pulsation Therapy



Sapphire Pulse™ features pulsation and static low air loss therapy modes for pressure injury prevention. Pulsating cycle times can be adjusted from 3 to 20 minutes. Sapphire Pulse's high-strength air cells rest on a foam safety pad for patient protection during transfers, and the top cover is vapor-permeable and machine-washable.

- Pulsation and Static therapy modes
- Cycle times adjustable from 3 to 20 min.
- Vapor-permeable, machine-washable top cover
- Fowler feature increases airflow when head of bed is elevated

Key Specs

MATTRESS

Capacity 600 lb./1,000 lb.
Width 35", 39", 42", or 48"
Length 80"

CONTROL UNIT

Zoned Low Air Loss Yes
Pulsating Low Air Loss 3-20 min.
Alarms Power Failure and Low Pressure
Pressure Sensor Technology Yes

FABRIC OPTIONS

Nylon 70-Denier–Polyurethane-Coated .. Standard
4-Way Stretch Optional
Non-Skid Polyurethane Bottom Cover ... Standard

FABRIC CHARACTERISTICS



Powered
Surface



Cool
Comfort



Low
Friction

sizewise.com 800-814-9389