

Sapphire Static LAL™

Low Air Loss Therapy



Key Specs

MATTRESS

Capacity 600 lbs.
Width 35"
Length 80"

CONTROL UNIT

Alarms Power Failure and Low Pressure
Pressure Sensor Technology Yes
Clinically Effective Low Air Loss Yes

FABRIC OPTIONS

Nylon 70-Denier–Polyurethane-Coated .. Standard
PVC Laminate Bottom Cover Standard

Sapphire Static LAL™ is a complete mattress replacement system with three zones for optimal low air loss therapy. The top cover is waterproof, vapor-permeable, and machine-washable. Diffuses 1,740 liters of air per minute for patient dryness and comfort.

- *Three zones for low air loss therapy*
- *Waterproof, vapor-permeable, machine-washable top cover*
- *1,740 liters of air per minute*

size wise.com 800-814-9389