

Shuttle™ A and B Series

Safe Patient Handling and Mobility Chair



Shuttle A (Standard)

Shuttle B (Bariatric)

Shuttle™ A and B Series wheelchairs transform from chair to stretcher, and are ideal for use with lateral transfer devices. Shuttle's leg supports elevate independently for enhanced patient comfort, and the chair features a stand-assist function and ambulation poles to aid caregivers with patient mobility.

- *Transforms from stretcher to chair to stand-assist device*
- *Ideal for use in lateral transfers*
- *Leg supports elevate independently*
- *Standard (Shuttle A) and Bariatric (Shuttle B) versions available*
- *Optional Power Drive enhances safe patient handling*

sizewise.com 800-814-9389


Sizewise®
Designed to heal



Key Specs

SHUTTLE A SERIES

Weight Capacity	650 lbs.
Seat Width	25"
Seat Depth	18"
Seat Height from Floor (low)	23.5"
Seat Height from Floor (high)	31.5"
Stretcher Height from Floor (low)	29.5"
Stretcher Height from Floor (high)	37.5"
Degree of Tilt	15°
Power Drive	Optional

SHUTTLE B SERIES

Weight Capacity	1,000 lbs.
Seat Width	31"
Seat Depth	18"
Seat Height from Floor (low)	23.5"
Seat Height from Floor (high)	31.5"
Stretcher Height from Floor (low)	29.5"
Stretcher Height from Floor (high)	37.5"
Degree of Tilt	15°
Power Drive	Optional

- 1 — Shuttle features three modes: chair, stretcher, and stand-assist
- 2 — Leg supports elevate independently for patient comfort
- 3 — Optional Power Drive enhances safe patient handling

sizewise.com 800-814-9389

