

# Sleep Lab Bed™

*Air-Adjustable Support Surface  
with Medical-Grade Top Cover*



*Shown with optional  
Adjustable Base™*

Sleep Lab Bed features durable materials and air-adjustable comfort to support the unique needs of a busy sleep center. The mattress offers 45 levels of soft/firm adjustability to help provide patients with optimal sleep conditions during their study, while the wipeable top cover and RF-welded seams help protect the integrity of the surface for a longer useful life. Different mattress sizes and frame options are available to provide sleep labs with greater flexibility in addressing a wide range of patients.



*Medical-grade design and construction  
with consumer-grade features and comfort*



*Performance fabrics, RF-welded seams  
and Core Shield™ improve durability and  
minimize risk of fluid ingress*



*CuPro™ copper- and gel-infused foam helps  
enhance microclimate by dissipating heat  
and moisture*

## MATTRESS

Max Patient Weight ..... 1,000 lb.  
Width x Length (Full XL)..... 54" x 80"  
Width x Length (Queen)..... 60" x 80"  
Height ..... 10.5"

Comfort Level Settings ..... 45  
Wired Remote ..... Yes  
Memory Settings..... 1  
Bases ..... See Reverse

## PERFORMANCE FOAM AND FABRICS

- Polycarbonate 4-way stretch top cover
- CuPro™ copper- and gel-infused foam
- Core Shield™ protective layer
- Non-skid polyurethane bottom cover



**LEARN MORE**

Ask your Agiliti representative or call 800-814-9389




## ADJUSTABLE BASE™



Our optional Adjustable Base is built with medical-grade, high-density polyethylene—making it extra durable and easy to clean.

- 850 lb. capacity
- Sizes: Full XL or Queen
- Wireless hand control
- Lockable 4" casters for stability and safety
- Optional siderails and headboard
- UL 962 compliant: household and commercial use



## OTHER BASE OPTION



### Mighty Rest®

- 1,000 lb. capacity
- Fully adjustable medical bed frame
- Sizes: Full XL or Queen
- IEC 60601 Certified