

# Stratus™ 7.1

Low Air Loss Alternation Therapy



## Key Specs

### MATTRESS

Capacity (Standard) .....	550 lbs.
Capacity (Bariatric) .....	650 lbs.
Width .....	35", 42"
Length .....	80"

### CONTROL UNIT

Liters of Air Per Minute .....	12
Cycle Times .....	10-25 min
Zoned Air Loss .....	Yes
Alarms .....	Low Pressure, Power/System Failure
Quick Inflator .....	Optional

### FABRIC OPTIONS

4-Way Stretch Top Cover .....	Standard
PVC Laminate Bottom Cover .....	Standard

Stratus 7.1™ offers clinically effective 1-in-3 alternation therapy for the prevention and treatment of pressure injuries. The mattress features 20 individual high-strength air cushions for pressure redistribution and blood flow stimulation, powered by an easy-to-use digital pump.

- *1-in-3 alternation therapy for the prevention and treatment of pressure injuries*
- *Variable cycle times: 10-25 min.*
- *Whisper-quiet digital pump*

size wise.com 800-814-9389