

Stratus™ 9.1

Low Air Loss Alternation and Pulsation Therapy



Stratus 9.1™ offers clinically effective 1-in-3 alternation therapy with simultaneous pulsation for the prevention and treatment of pressure injuries. The mattress features 20 individual high-strength air cells for pressure redistribution and blood flow stimulation, powered by an easy-to-use blower.

- *Three therapy modes: 1-in-3 alternation with pulsation, 1-in-3 alternation, and pulsation*
- *Indicated for the prevention and treatment of pressure injuries*
- *Microprocessor-controlled blower with variable cycle times: 3-20 min.*
- *Includes static, fowler, and autofirm features*

Key Specs

MATTRESS

Capacity (Standard) 600 lbs.
Capacity (Bariatric) 1,000 lbs.
Width 35", 39", 42", 54" or 60"
Length 80"
Standard and Bariatric Sizes Available

CONTROL UNIT

Liters of Air Per Minute 1,740
Zoned Air Loss Yes
Alarms Power Failure and Low Pressure

FABRIC OPTIONS

Nylon 70 Polyurethane-Coated Top Cover Standard
4-Way Stretch Top Cover Optional
Non-Skid Polyurethane Bottom Cover Standard

sizewise.com 800-814-9389


Sizewise®
Designed to heal