

Tom 2

Pediatric Crib by LINET®




Shown with Pediatric Pulse™ LAL support surface. Also available with foam.

Tom 2 Pediatric Crib offers innovative features that improve the overall patient experience. Designed with a focus on safety, Tom 2 supports fall prevention protocols, improves accessibility to the patient, and provides a safe work environment for caregivers. Telescopic side rails with five height positions provide a customizable solution for safe and easy access to the child.

- Easy, one-button height, backrest, and trend/reverse trend adjustments
- Telescopic side rails with five positions provide easy patient access
- Built-in scale and multi-zone bed exit alarm
- Modern design includes transparent bed ends, giving children a better view of their surroundings

sizeWise.com 800-814-9389


SizeWise®
Designed to heal



5

Sizewise® Support Surface Options + LINET® Bed Frame Technology = **Sizewise®**
Tom 2 by LINET®

Key Specs

BED

Low Position – Deck to Floor	23"
High Position – Deck to Floor	33"
Deck Width	28"
Deck Length	56"
Safe Working Load	165 lbs.

SIDE RAILS

Type	Telescopic
Positions	5
High Position	32"
Middle Position	16"
Low Position	0"

CASTER AND BRAKING SYSTEM

Central Caster Control	Standard
----------------------------------	----------

BED ANGLES

Backrest/Head Section	0-40°
Trend	15°
Reverse Trend	-15°

OPTIONS AND ACCESSORIES

Bed Exit Alarm	Standard
One-Touch Bed Positioning	Standard
Scale	Standard

MATTRESS OPTIONS

Tom 2 pairs clinically with Pediatric Pulse™ and standard foam support surfaces.

- 1 — Telescopic side rails stop at five positions for enhanced access
- 2 — Autoregression feature ensures a constant safety gap between backrest and head of crib
- 3 — Transparent bed ends give children a better view of their surroundings
- 4 — Available with Pediatric Pulse™ LAL pulsation support surface
- 5 — A unique cross-manufacturing agreement provides the best of both worlds: Sizewise® surface options and LINET® bed frames

sizewise.com 800-814-9389

