

# Victory 500™

Foam Support Surface



Victory 500™ foam support surface provides pressure redistribution and shear reduction for optimal patient comfort. High-resiliency viscoelastic foot end is softer than the rest of the surface to protect vulnerable heels. Full-length base foam helps prevent bottoming out. Optional Mattress Asset Protection (MAP) Liner™ safeguards internal components and reduces shear.

- *Low-shear, antimicrobial, fluid-resistant top cover*
- *Full-length base foam prevents bottoming out*
- *Optional RF-welded seams for max infection control*
- *Made-to-order item*

## Key Specs

Capacity .....	450 lbs./650 lbs.
Width .....	.35", 39", or 42"
Length .....	.76", 80", or 84"
Height.....	7"

## CONSTRUCTION

Foam Layers.....	3
MAP Liner™ .....	Optional
Fire Barrier .....	Standard
Defined Perimeter.....	Optional

## FABRIC OPTIONS

Envy® 4-Way Stretch Top Cover .....	Standard
Envy® Non-Skid Polyurethane .....	Standard
RF-Welded Seams .....	Optional
Polycarbonate 4-Way Stretch Top Cover . . .	Optional

sizewise.com 800-814-9389

