

Victory 900™

Pressure Redistribution Support Surface



Key Specs

Capacity	450 lbs./650 lbs.
Width	35", 39", or 42"
Length	76", 80", or 84"
Height	7"

CONSTRUCTION

Foam Layers	5
MAP™ Liner	Optional
Fire Barrier	Standard
Defined Perimeter	Optional

FABRIC OPTIONS

Envy® 4-Way Stretch Top Cover	Standard
Envy® Non-Skid Polyurethane	Standard
RF-Welded Seams	Optional
Polycarbonate 4-Way Stretch Top Cover	Optional

Victory 900™ features four layers of impression load deflection (ILD) foam, housed in a 1/2" viscoelastic memory-style foam for optimal patient comfort. The viscoelastic foot end is softer than the rest of the surface to protect vulnerable heels, and a concealed zipper prevents skin tears and protects inner foam.

- *Pressure redistribution surface*
- *Heel area softer than rest of surface*
- *Optional RF-welded seams for max infection control*
- *Made-to-order item*

sizeWise.com 800-814-9389

