



Efficiently Order Medical Equipment Directly from Your EMR System

Save time for your clinical staff. Enable them to place rental orders for specialty beds, surfaces and other medical equipment directly from your Epic or Cerner EMR system. The EMR-to-Agiliti ordering integration provides familiar workflows that remove ordering bottlenecks and improve productivity.

Automate the ordering process

Order the exact equipment patients need and attach relevant patient information.

- Utilize just-in-time orders for improved time to therapy
- Reduce opportunities for ordering errors

Trigger automated pickup

Configure to automatically send pickup requests to Agiliti upon patient discharge.

- Reduce nurse burdens with fewer administration tasks

Technical requirements and effort

- Integrates with EPIC and Cerner EMRs
- Facility/IDN IT team allows orders to send via HL7 v.2
 - Connect through one of the Agiliti standard patterns:
 - TCP/IP
 - SFTP
 - REST
 - FHIR
- Customized but on average between 20-60 hours of integration work from facility IT or EMR admin to build an equipment catalog, order panel and Agiliti interface.
 - Typically performed over 6-8 weeks
 - Proprietary EMR-to-Agiliti build guide reference provided
 - Agiliti offers outsourced resources for integration work upon request

DRIVE EFFICIENCIES FOR YOUR CLINICAL TEAM



REDUCE
TRAINING NEEDS



IMPROVE
CAREGIVER PRODUCTIVITY



REMOVE
ORDERING BOTTLENECKS