

# QUICK-SHIP FOAM SUPPORT SURFACES

Agiliti quick-ship options provide facilities a fast, cost-effective way to replace damaged foam support surfaces. Products can be delivered in limited quantities within five business days. Contact your Agiliti representative for production timeframe on larger quantities.



NP Trio™



NP Adjust™



NPT4™

<b>Product Description</b>	3-layer, 3-zone non-powered foam support surface	Vented, self-adjusting foam support surface	Vented, self-adjusting foam support surface with optional pump
<b>Standard Size</b>	35"W x 82"L	35"W x 82"L	35"W x 82"L
<b>Max Patient Weight</b>	500 lb.	500 lb.	500 lb.
<b>Warranty: Top/Bottom Cover*</b>	4 Years	4 Years	2 Years
<b>Warranty: RF-Welded Seams*</b>	4 Years	4 Years	5 Years
<b>Product Webpage</b>	<a href="https://agilitihealth.com/product/np-trio">agilitihealth.com/product/np-trio</a>	<a href="https://agilitihealth.com/product/np-adjust">agilitihealth.com/product/np-adjust</a>	<a href="https://agilitihealth.com/product/npt4">agilitihealth.com/product/npt4</a>

These Agiliti-manufactured support surfaces feature construction and design details that help serve patients in a busy hospital environment:



**4-way stretch top cover** fabric is durable, breathable, waterproof and highly chemically resistant



**CoreShield™ liner** provides an additional layer of waterproof protection to the foam core (MAP Liner™ on NPT4)



**RF-welded seams** protect against fluid ingress, extending the useful life of the product



**CuPro™ copper- and gel-infused foam** helps dissipate heat from the patient/surface interface

## Learn more:

Contact your Agiliti sales representative or call 800-814-9389