

Approved Cleaners (or equivalent):	
Product Name	Manufacturer
Clorox Healthcare® Bleach Germicidal Wipes	The Clorox Company
Clorox® Dispatch Bleach Germicidal Cleaner	The Clorox Company
Expose® II 256	Diversey, Inc.
Sani-Cloth® Bleach Germicidal Disposable Wipes	Nice-Pak/PDI, Inc.
MICRO-KILL™ Bleach Germicidal Wipes	Medline Industries, Inc.
MICRO-KILL™+ Disinfectant Wipes	Medline Industries, Inc.
Oxivir® 1 Wipes	Diversey, Inc.
Oxivir® Five 16 Concentrate (US) (properly diluted)	Diversey, Inc.
Oxivir® TB Wipes Based	Diversey, Inc.
Oxycide™ Daily Disinfectant Cleaner	Ecolab Inc.
PeridoxRTU® Sporocidal Disinfectant and Cleaner	Contec, Inc.
Rubbing Alcohol (99%)	Not applicable
Sani-Cloth® Plus Germicidal Disposable Cloth	Nice-Pak/PDI, Inc.
Super Sani-Cloth® Germicidal Wipe	Nice-Pak/PDI, Inc.
Virex® II 256 One-Step Disinfectant Cleaner and Deoderant	Diversey, Inc.
Unsupported Cleaners*:	
Product Name	Manufacturer
Clorox Commercial Solutions® Clorox® Disinfecting Wipes	Clorox Professional Products Company
Lysol® Brand II Disinfecting Wipes	Reckitt Benckiser Inc.
MICRO-KILL® One Germicidal Alcohol Wipes	Medline Industries, Inc.
Quat Plus TB	Medline Industries, Inc.
Sani-Cloth® AF3 Germicidal Disposable Wipes	Nice-Pak/PDI, Inc.
Virex® TB	Diversey, Inc.

* Internal test results have shown that the cleaners listed will reduce some plastic component life with repeated use. Products sanitized with these cleaners are not covered by manufacturer warranty.

Warning

Some cleaning agents will reduce the life of some plastic components with repeated use. Before using the system check all plastic components for those that may become brittle and may break under stress. This could result in a delay in therapy. See manufacturer list of approved cleaners for more information.

For additional product support information, visit:
agilitihealth.com/clinical-support/patient-care-plus-resources