


PRODUCT SPOTLIGHT: ESSENTIA™ MULTI-ACUITY BED FRAME

Unique Product Features Enhance Clinical Versatility

Essentia™ provides hospitals one versatile platform to serve patients across all acuity levels (Med-Surg to ICU), with therapeutic support surface options that can be tailored per patient risk profile. New, innovative features – including the first-of-its-kind SlideRail™ – help lower caregiver risk and support patient mobility.



AGILITI EXCLUSIVE

SlideRail™

- Industry-first design easily slides toward the foot of the bed
- 39" unobstructed ingress/egress and multiple handholds make it easier for patients to safely exit or enter the bed



AGILITI EXCLUSIVE

RevTrend 20™

with Weight-bearing Footboard

- 20-degree Reverse Trendelenburg position
- Weight-bearing footboard allows progressive strengthening exercises



AGILITI EXCLUSIVE

MotionStop 360™

- 360-degree sensors around the base stop bed from lowering if an obstruction is detected



Siderail Height

- Accommodates therapeutic support surfaces up to 10.5" thick for greater immersion
- Design aligns with patient safety guidelines (60601-2-52)



PowerBrake™ with AutoLock™

- Brakes automatically lock after 60 seconds while bed is plugged in
- Brake/Neutral/Steer can be set with the press of a button at three locations around the bed



SafeGlance™

- Green lights, visible from a distance, indicate whether bed is in its safest state
- Shows status of bed exit alarm, protocol timer alarm, brakes and bed low position

Learn More

Ask your Agiliti representative or call 800-814-9389