



## Equipment Location and Status Tracking

inCare™ Equipment Tracking Software is the key proprietary technological component of the Agiliti Onsite Management Program.

### How It Works

Each piece of equipment is assigned a unique barcode, allowing Agiliti onsite staff to scan and update equipment status information through a handheld scanner. The ability to track and monitor equipment empowers Agiliti staff to identify utilization trends, locate missing equipment, minimize cross-contamination risk, manage safety inspections and maintenance, and much more.

### Customizable Reporting

Customizable reporting through inCare allows Agiliti onsite staff to track equipment productivity, resulting in an accurate portrayal of utilization and elimination of excess equipment expenses. Designed with HIPAA and The Joint Commission guidelines in mind, this system helps maintain regulatory compliance in your facility.

### Agiliti Will Provide

- inCare software application
- SQL database
- Barcode scanner
- ADT Receiving Engine



**THE RIGHT  
EQUIPMENT,  
in the right place,  
at the right time.**



## Implementation Components:

### Hospital Server

(physical or virtual)

- Windows Server OS
- Microsoft SQL Server

(Agiliti provided)

- HL7 receiver service
- HL7 interpreter service

### Supported Operating Systems

- Windows 10
- Windows Server 2012 or higher

### Database Server

- Microsoft SQL Server 2012 or higher
- Core i5 or higher
- 4 GB ram
- 120 GB hard disk
- 2 GB allocated space for database



### Workstation PC, Monitor, and Printer

(hospital or Agiliti provided)

### Workstation PC Requirements

- Windows 10
- Core i5 or higher
- 4 GB recommended
- 80 GB hard disk
- USB port



Network Connection

### Handheld Barcode Scanner

(Agiliti provided)

- inCare mobile application



USB Connection

Honeywell CT40

## Overall Requirements:

### Required Programs:

- Microsoft Excel 2010 (or higher)
- Adobe Reader/PDF creator

### Programs Provided by Agiliti:

- MDLink Receiver
- Agiliti Custom HL7 Processor (local or server)
- Bomgar Web Meeting

### Required Ports

- Port 992 - Outbound only to rocky.agilitihealth.com
- Port 22 - Inbound and Outbound to rocky.agilitihealth.com

### Admits, Discharges and Transfers (ADT)

Facilities may enter information by hand, but most choose to use a direct method of sending ADT messages from the hospital to inCare. ADT messages are sent directly from the EMR to inCare in real-time with no action required on the part of the onsite program staff.

Contact your Agiliti Representative to learn more about inCare™ Equipment Tracking System.

### About Agiliti

Agiliti helps hospitals reclaim millions of dollars in wasted medical equipment. The company places equipment management experts and clinical engineers in more than 7,000 hospitals to ensure that the right equipment is in the right place at the right time. Onsite services free nurses and other clinicians to spend more time caring for patients and less time searching for the equipment they need. Agiliti also provides responsive rental services for a full range of equipment, including respiratory systems, infusion therapy devices, specialty beds and surgical lasers. A nationwide network allows the company to serve all health systems and hospitals across the U.S.

For more information, visit [agilitihealth.com](http://agilitihealth.com).

### CORPORATE OFFICE

6625 West 78th Street. Suite 300. Minneapolis, MN. 55439

Toll Free: 800.847.7368 Phone: 952.893.3200 Fax: 952.893.0704 [www.agilitihealth.com](http://www.agilitihealth.com)

