










PRODUCT	DESCRIPTION
 <p>BIOJET™ MRI ULTRASOUND FUSION GUIDED PROSTATE BIOPSY</p>	<p>Enhanced lesion targeting by combining MRI and real-time ultrasound for prostate biopsies</p> <ul style="list-style-type: none"> • Support for both transrectal and transperineal biopsy procedures • Studies support the identification of 50% more prostate cancer findings • Available for use in office, hospital or Ambulatory Surgery Center (ASC)
 <p>BK FLEX FOCUS 800 and 5000 ULTRASOUND</p>	<p>Robotic, laparoscopic and hand-assist partial nephrectomy or prostate biopsies</p> <ul style="list-style-type: none"> • Real-time image for tumor margin identification and complete removal • Excellent color Doppler with superb spatial resolution and sensitivity
 <p>Aloka® ProSound Alpha 6SX Ultrasound</p>	<p>Robotic, laparoscopic and hand-assist partial nephrectomy or prostate biopsies</p> <ul style="list-style-type: none"> • Real-time image for tumor margin identification and complete removal • Doppler capability to visualize and account for blood flow • High resolution digital imaging
 <p>Storz Modulith® SLX and SLX-F2 MOBILE ESWL</p>	<p>Extracorporeal shockwave lithotripsy to treat kidney and ureteral calculi</p> <ul style="list-style-type: none"> • Patented cylinder source provides unsurpassed fragmentation efficiency with low side effects • Focal depth 165-180mm • 300-500 lbs. patient weight capacity • Endo-urology table with Trendelenburg
 <p>Dornier® Compact Delta 2 Mobile ESWL</p>	<p>Extracorporeal shockwave lithotripsy to treat kidney and ureteral calculi</p> <ul style="list-style-type: none"> • Integrated shock source and C-arm to make setup fast and reliable • Isocentric movements ensure precise targeting • Unique, motorized therapy head allows over and under table treatments
 <p>BOSTON SCIENTIFIC/AMS Greenlight XPS™ Laser</p>	<p>Photo-vaporization of the prostate</p> <ul style="list-style-type: none"> • 180 watts of 532nm laser energy with liquid cooled MoXy™ fiber • Minimally invasive outpatient procedure • TruCoag to cauterize vessels
 <p>EndocarE® CryoCarE® CS CryoTherapy System</p>	<p>Cryosurgical ablation of prostate and kidney cancer</p> <ul style="list-style-type: none"> • Percutaneous, laparoscopic or open renal cryoablation • Variable probe technology to precisely treat cancerous lesions • Planning software to treat various cancers



PRODUCT	DESCRIPTION
 <p>Galil Visual-ICE® Precise and SeedNet CryoTherapy System</p>	<p>Cryosurgical ablation of prostate and kidney cancer</p> <ul style="list-style-type: none"> • Percutaneous, laparoscopic or open renal cryoablation • Innovative needle design • Planning and simulation software to treat various cancers
 <p>Boston Scientific/AMS Swiss LithoClast® SELECT Intracorporeal Lithotripter</p>	<p>Bladder, kidney, ureteral calculi removal</p> <ul style="list-style-type: none"> • Useful in fragmentation or evacuation of calculi while simultaneously using ultrasonic and pneumatic lithotripter • New Vario ultrasound handpieces are 28% lighter weight designed for more comfort and convenient handling • New suction adapter angle designed to improve clearance
 <p>Lumenis® VersaPulse PowerSuite (20, 60, 100 Watt)</p>	<p>Ureteral, kidney and bladder stone lithotripsy</p> <ul style="list-style-type: none"> • Fragmentation of any size stones in urinary tract • Efficient vaporization or enucleation of the prostate • Vaporization and coagulation of soft tissue, strictures and bladder tumors
 <p>Revolix® 200 Watt Thulium Laser</p>	<p>Vaporization, excision, and coagulation of soft tissue, prostate strictures, bladder tumors and partial nephrectomy</p> <ul style="list-style-type: none"> • Excellent visualization and speed for the enucleation or ablation of prostate • Superior vaporization, cutting and hemostasis of prostate tissue • Eliminates the risk of unintentional damage to adjacent tissue
 <p>Richard Wolf® Piranha Morcellator</p>	<p>Removes prostatic tissue with speed and efficiency after laser enucleation of the prostate</p> <ul style="list-style-type: none"> • Quick and safe removal of large adenomas • Great visualization which significantly reduces operating time with 50% inlet and 50% outlet

About Agiliti

Agiliti helps hospitals reclaim millions of dollars in wasted medical equipment. The company places equipment management experts and clinical engineers in more than 7,000 hospitals to ensure that the right equipment is in the right place at the right time. Onsite services free nurses and other clinicians to spend more time caring for patients and less time searching for the equipment they need. Agiliti also provides responsive rental services for a full range of equipment, including respiratory systems, infusion therapy devices, specialty beds and surgical lasers. A nationwide network allows the company to serve all health systems and hospitals across the U.S.

For more information, visit agilitihealth.com.

CORPORATE OFFICE

6625 West 78th Street. Suite 300. Minneapolis, MN. 55439

Toll Free: 800.847.7368 Phone: 952.893.3200 Fax: 952.893.0704 www.agilitihealth.com

