



- ▶ Ultra-low deck height minimizes fall-related injuries
- ▶ Bed exit alarm and innovative triple-guard anti-entrapment system provide optimum patient safety
- ▶ Electronic central brake and steer for safe and easy patient transport
- ▶ Antimicrobial surfaces support infection control initiatives

Harmony® Low Bed: A pitch-perfect blend of patient care and innovation

The OnCare™ Harmony® Low Bed rises above typical acute care beds with a long list of innovative standard features to help caregivers achieve optimum patient safety. A minimum deck height below 10", bed exit alarm and triple-guard anti-entrapment system minimize patient injuries. In addition, easy-to-use controls, such as the one-touch electronic central brake and steer, let caregivers focus on patient care. These features give you a smarter, better-performing option for acute care: the Harmony Low Bed.

Meeting your standard of care

The Harmony Low Bed was built on insights gained from working alongside caregivers, in thousands of facilities, every day for more than 70 years.

Acquisition made affordable

OnCare works with you to make equipment acquisition affordable through rental, purchase or lease.



Harmony® Low Bed Features and Specifications



Standard Features

- Integrated bed exit alarm
- Integrated patient scale
- Nurse call integration
- Electronic central brake and steer
- Triple-guard anti-entrapment system
- Manual and electronic CPR release
- Transportability at any height
- Antimicrobial plastic surfaces
- Under bed light
- Integrated USB port
- Fully electric articulation
- Integrated IV pole
- Integrated pump holder
- Bed length extension
- Auto contour
- Auxiliary power outlet
- Patient control hand pendant

Optional

- Trapeze/Patient helper



Electronic central brake and steer



Nurse call and USB port



Anti-entrapment sensors



Backlit, intuitive controls

Specifications

Deck Height:	
Low	9.875"
High	31.5"
Overall Length	91.75"
Overall Width	40.5"
Sleep Deck	35" x 80" (expandable to 88")
Weight Capacity	400 lbs.
Head Section Adjustment	0°-70°
Knee Gatch Adjustment	0°-35°
Trendelenburg and Reverse Trendelenburg	14°
Caster Diameter	5"

Electronics

Volts	120 VAC
Ampere Rating	5.0 Amps @ 120 VAC 60Hz or 600 VA (Max Power Limit)
Current Leakage	<10mA (RMS)

Note: Specifications subject to change without notice.

Visit www.oncaremedical.com to learn more.

OnCare™ Harmony® is a registered trademark of Agiliti Health, Inc.

CORPORATE OFFICE

6625 West 78th Street, Suite 300, Minneapolis, MN, 55439

Toll Free: 800.847.7368 Phone: 952.893.3200 Fax: 952.893.0704 www.agilitihealth.com

©2018 Agiliti Health, Inc. HARMONY1218

